

UDITracker® OR RFID Storage

A SEAMLESS STORAGE AND TRACKING SOLUTION

Our RFID storage fully integrates into UDITracker® to completely automate your chain of custody and device management. The Swipe, Grab and Go!™ access, ensures that all transactions are automatically captured. This minimizes data entry, reducing retrieval time and increasing staff accountability, so you can focus less on paperwork and more on patients.

- No data entry
- Improve regulatory compliance
- Control cost and waste

Call us today to get started!

866.803.3720 | www.ChampionHT.com

Refrigerator
Model TS029

“We love them.
Super easy to use
and really makes
tracking so much
more efficient.”

Debbie Ramos
Hendrick Medical Center
Abilene, TX

CHAMPION
HEALTHCARE TECHNOLOGIES

Freedom to do more good.

REFRIGERATOR MODEL TS029

Swipe, Grab & Go!™

Using secure access, each inventory transaction is automatically captured and processed in real-time.

1. Place internet enabled RFID cabinet at the point of use
2. Pre-authorized users access the cabinet with personalized badges
3. Users add or remove RFID tagged inventory from the cabinet



Volume: 5.5 ft³

External Dimensions: 42" x 24" x 26"

Internal Dimensions: 28" x 21" x 18"

Additional options

- Open shelves
- Slide trays with dividers

Specific applications

- Orthopedic implants
- Cardiovascular implants
- Neurosurgery implants
- Tissue grafts

REFRIGERATOR MODEL TS029

Physical Characteristics	Imperial	Metric
External Dimensions	42h x 24w x 26d in	1054h x 610w x 660d mm
Internal Dimensions	28h x 21w x 18d in	711h x 533w x 457d mm
External Space	42h x 24w x 47d in Place 3in from back wall and 2in from front wall	1054h x 610w x 1194d mm Place 76mm from back wall and 534mm from front wall
Volume	14.9 ft ³ (useable)	156 L
Weight	165 lbs	75 kg
Control Range	36° F to 46° F	2° C to 8° C
Security Device	Electromagnetic	
Communications	USA	
RFID Frequencies	902 to 928 MHz	
Communications	Ethernet based LAN DSL optional Cellular optional	
Access Control	EM4102/EM4200 Proximity Reader	
Power		
Input Rating	120VAC, 60Hz, 8A	
Output Rating	120VAC, 60Hz, 5A* max (*only if connected device has its own power protection)	
Power Consumption	160W max	
Battery	Two - 12VDC, 1.3Ah, Sealed lead acid	
Refrigerant	R - 134a	
Cooling Method	Forced Air Circulation	
Environmental		
Environmental Protection	IPX0	



Cabinets provided by Champion Healthcare Technologies in partnership with Terso Solutions, inc.

Champion Healthcare Technologies is committed to providing reliable RFID solutions that are safe for those who use them. Therefore, all of our products are certified by the FCC for EMC emissions and have received UL Certifications.