

# UDITracker® OR RFID Storage

## A SEAMLESS STORAGE AND TRACKING SOLUTION

Our RFID storage fully integrates into UDITracker® to completely automate your chain of custody and device management. The Swipe, Grab and Go!™ access, ensures that all transactions are automatically captured. This minimizes data entry, reducing retrieval time and increasing staff accountability, so you can focus less on paperwork and more on patients.

- No data entry
- Improve regulatory compliance
- Control cost and waste

**Call us today** to get started!

866.803.3720 | [www.ChampionHT.com](http://www.ChampionHT.com)

Ultra Low  
Temp Freezer  
Model TS057

“Before our UDITracker® Freezer, we were writing off up to \$200,000 per year in missing tissue grafts. My Director is so pleased with the results that we will be adding additional units when our new OR is completed.”

**Argenis Mendoza**  
Perioperative Information Systems Analyst  
Reston Hospital Center

**CHAMPION**  
HEALTHCARE TECHNOLOGIES

Freedom to do more good.

# ULTRA LOW TEMP FREEZER MODEL TS057

## Swipe, Grab & Go!™

Using secure access, each inventory transaction is automatically captured and processed in real-time.

1. Place internet enabled RFID cabinet at the point of use
2. Pre-authorized users access the cabinet with personalized badges
3. Users add or remove RFID tagged inventory from the cabinet

Swipe,



Grab & Go!™



**Volume:** 14.9 ft<sup>3</sup>

**External Dimensions:** 78" x 27" x 37.6"

**Internal Dimensions:** 28" x 21" x 18"

### Additional options

- Open shelves
- Slide trays with dividers

# ULTRA LOW TEMP FREEZER MODEL TS057

Physical Characteristics	Imperial	Metric
External Dimensions	78h x 27w x 37.6d in	198.1h x 68.6w x 95.5d cm
Internal Dimensions	28h x 21w x 18d in	130.1h x 45.1w x 71.86d cm
External Space	Place 3in from back wall and 2in from front wall	Place 76mm from back wall and 534 mm from front wall
Volume	14.9 ft <sup>3</sup> (useable)	421 L
Weight	165 lbs	75 kg
Control Range	-58° F to -122° F	-50° C to -86° C
Security Device	Electromagnetic	
Communications	USA	
RFID Frequencies	902 to 928 MHz	
Communications	Ethernet based LAN DSL optional Cellular optional	
Access Control	EM4102/EM4200 Proximity Reader	
Power		
Input Rating	208VAC, 60Hz, 8A 1x NEMA 5-20/1x NEMA -15 (requires single dedicated circuit)	
Output Rating	208-230VAC	
Power Consumption	1.9KW max	
Battery	Two - 12VDC, 1.3Ah, Sealed lead acid	
Refrigerant	R - 134a	
Cooling Method	Forced Air Circulation	
Backup System	LN2 or CO2 Back-up System	
Environmental		
Environmental Protection	IPXO	



Cabinets provided by Champion Healthcare Technologies in partnership with Terso Solutions, inc.

Champion Healthcare Technologies is committed to providing reliable RFID solutions that are safe for those who use them. Therefore, all of our products are certified by the FCC for EMC emissions and have received UL Certifications.